

## Detailed information

|                   |   |
|-------------------|---|
| Material          | Synthetic Nitrile latex.  |
| Type              | Non-Sterile Powder-Free.<br>Ambidextrous; Finger Tip Textured; Beaded Cuff; White or Coloured (Blue, Light Blue, ...)   |
| Quality Standards | Conforms to ASTM D6319<br>Manufactured under ISO9001: 2008, ISO 13485:2003. ISO 22000:2005 Quality Management System.<br>Manufactured from 100% nitrile (Acrylonitrile-Butadiene) |
| Gloves Size       | Extra-small, Small, Medium, Large, Extra-large.<br>Marked in the check box on the shipping carton with black ink.   |
| Storage           | Store in a dry and cool place, the temperature not higher than 38°C.  |
| Shelf-life        | 3 years from the date of manufacturing.   |

 DISHANG

# DISPOSABLE NITRILE GLOVES

ANTI-CORROSION / OIL-PROOF / MORE COMFORT / MORE SECURITY



HOME



MEDICAL



EXPERIMENT



FACTORY





## PARAMETER INFORMATION

|               |  |
|---------------|--|
| <b>NAME</b>   | DISPOSABLE NITRILE GLOVES              |
| <b>COLOUR</b> | BLUE (COMMON FOR LEFT AND RIGHT HANDS) |
| <b>APPLY</b>  | MEDICAL EXAMINATION GLOVES (CLASS I)   |

**TYPE** NO POWDER / DISPOSABLE

**SPECIFICATIONS** 100PCS/BOX

**TEXTURE OF MATERIAL** NITRILE

\*NOTE: IF THERE IS ANY ERROR IN THE INSPECTION DATA, PLEASE REFER TO THE REAL OBJECT

**SIZE** S: LONG 23.5CM WIDE 8CM M: LONG 24CM WIDE 9CM  
L: LONG 23.5CM WIDE 10CM XL: LONG 24CM WIDE 11CM

## HOW TO SELECT SIZE

S SIZE

PALM WIDTH 7.1-8CM

M SIZE

PALM WIDTH 8.1-9CM

L SIZE

PALM WIDTH 9.1-10CM

XL SIZE

PALM WIDTH 10.1-11CM

MOST WOMEN USE S, AND THOSE WITH LARGER HANDS USE M  
MOST MEN USE SIZE M, AND THOSE WITH LARGER HANDS USE SIZE L

### REMINDER

DUE TO DIFFERENT MEASURING METHODS OR TOOLS,  
NITRILE GLOVES AND LATEX PVC GLOVES FOR THE DIFFERENCE  
OF MATERIALS, THE SIZE ERROR OF ABOUT 2-5MM IS ALLOWED.







## WEARING COMMON SENSE



1. PLEASE TRIM YOUR NAILS BEFORE WEARING THEM. IF THE NAILS ARE TOO LONG OR SHARP, THE GLOVES WILL BE DAMAGED

2. WHEN WEARING, PLEASE USE THE BELLY PART OF YOUR FINGERS TO AVOID THE GLOVES SLIPPING

3. WHEN TAKING OFF GLOVES, FIRST TURN OVER THE GLOVES AT THE WRIST, AND THEN TAKE OFF THE GLOVES AT THE FINGERS

# VDFGH

---



OIL-PROOF &  
WATER-PROOF



RESISTANCE TO  
HYDROCHLORIC ACID



EXCELLENT  
ELASTICITY



FLEXIBLE &  
COMFORTABLE



ENVIRONMENTAL &  
NON-TOXIC



POWDER-FREE &  
HYPOALLERGENIC

**Hongray®**

DISPOSABLE VINYL/NITRILE BLEND EXAMINATION

**GLOVES**



NON-STERILE  
POWDER-FREE  
AMBIDEXTROUS  
SINGLE USE  
LATEX-FREE

Easy Donning Technology

XS  S  M  L  XL

CE

**100**  
GLOVES



DISPOSABLE VINYL/NITRILE BLEND EXAMINATION  
**GLOVES**

NON-  
STERILE

POWDER-  
FREE

AMBIDE-  
XTROUS

SINGLE  
USE

LATEX-  
FREE

Easy Donning Technology



Vinyl/Nitrile Blend gloves containing no natural latex protein, causing no allergies, and offering affordable and effective safety and hygiene in various work environments for medical, industrial and food handling applications.

**Recommed product**



[Cleaning Towel manufacturer](#)



[Face towel](#)



[Cleaning Towel manufacturer](#)



[Microfiber beach towel](#)

If you are interested in our products, please feel free to contact us!!

Remember to look at our other [microfiber towel](#) or [beach towel](#). Thank you! ^^