

## Informazioni dettagliate

|                       |                                                                                                                                                                                       |
|-----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Materiale             | Lattice nitrile sintetico.                                                                                                                                                            |
| genere                | Senza polvere non sterile.<br>Ambidestro; Punta delle dita strutturata; Polsino in rilievo; Bianco o colorato (blu, azzurro, ...)                                                     |
| Standard di qualità   | Conforme a ASTM D6319<br>Prodotto secondo ISO9001: 2008, ISO 13485: 2003. ISO 22000: 2005 Sistema di gestione della qualità.<br>Prodotto al 100% in nitrile (acrilonitrile-butadiene) |
| Dimensione dei guanti | Extra-piccolo, Piccolo, Medio, Grande, Extra-grande.<br>Contrassegnato nella casella di controllo sulla confezione di spedizione con inchiostro nero.                                 |
| Conservazione         | Conservare in luogo fresco e asciutto, con temperatura non superiore a 38 ° C.                                                                                                        |
| Shelf-Fife            | 3 anni dalla data di produzione.                                                                                                                                                      |



# DISPOSABLE NITRILE GLOVES

ANTI-CORROSION / OIL-PROOF / MORE COMFORT / MORE SECURITY



HOME



MEDICAL



EXPERIMENT



FACTORY





## PARAMETER INFORMATION

|               |                                        |
|---------------|----------------------------------------|
| <b>NAME</b>   | DISPOSABLE NITRILE GLOVES              |
| <b>COLOUR</b> | BLUE (COMMON FOR LEFT AND RIGHT HANDS) |
| <b>APPLY</b>  | MEDICAL EXAMINATION GLOVES (CLASS I)   |

**TYPE** NO POWDER / DISPOSABLE

**SPECIFICATIONS** 100PCS/BOX

**TEXTURE OF MATERIAL** NITRILE

\*NOTE: IF THERE IS ANY ERROR IN THE INSPECTION DATA, PLEASE REFER TO THE REAL OBJECT

**SIZE** S: LONG 23.5CM WIDE 8CM M: LONG 24CM WIDE 9CM  
L: LONG 23.5CM WIDE 10CM XL: LONG 24CM WIDE 11CM

## HOW TO SELECT SIZE

S SIZE

PALM WIDTH 7.1-8CM

M SIZE

PALM WIDTH 8.1-9CM

L SIZE

PALM WIDTH 9.1-10CM

XL SIZE

PALM WIDTH 10.1-11CM

MOST WOMEN USE S, AND THOSE WITH LARGER HANDS USE M  
MOST MEN USE SIZE M, AND THOSE WITH LARGER HANDS USE SIZE L

## REMINDER

DUE TO DIFFERENT MEASURING METHODS OR TOOLS,  
NITRILE GLOVES AND LATEX PVC GLOVES FOR THE DIFFERENCE  
OF MATERIALS, THE SIZE ERROR OF ABOUT 2-5MM IS ALLOWED.







## WEARING COMMON SENSE



1.PLEASE TRIM YOUR NAILS BEFORE WEARING THEM. IF THE NAILS ARE TOO LONG OR SHARP, THE GLOVES WILL BE DAMAGED

2.WHEN WEARING, PLEASE USE THE BELLY PART OF YOUR FINGERS TO AVOID THE GLOVES SLIPPING

3.WHEN TAKING OFF GLOVES, FIRST TURN OVER THE GLOVES AT THE WRIST, AND THEN TAKE OFF THE GLOVES AT THE FINGERS

# VDFGH

---



OIL-PROOF &  
WATER-PROOF



RESISTANCE TO  
HYDROCHLORIC ACID



EXCELLENT  
ELASTICITY



FLEXIBLE &  
COMFORTABLE



ENVIRONMENTAL &  
NON-TOXIC



POWDER-FREE &  
HYPOALLERGENIC

**Hongray®**

DISPOSABLE VINYL/NITRILE BLEND EXAMINATION

**GLOVES**



Easy Donning Technology

NON-STERILE  
POWDER-FREE  
AMBIDEXTROUS  
SINGLE USE  
LATEX - FREE

XS ☐ S ☐ M ☐ L ☒ XL ☐

CE

**100**  
GLOVES



DISPOSABLE VINYL/NITRILE BLEND EXAMINATION  
**GLOVES**

NON-  
STERILE

POWDER-  
FREE

AMBIDE-  
XTROUS

SINGLE  
USE

LATEX-  
FREE

Easy Donning Technology



Guanti in vinile / nitrile che non contengono proteine del lattice naturale, non provocano allergie e offrono sicurezza ed igiene convenienti ed efficaci in vari ambienti di lavoro per applicazioni mediche, industriali e di manipolazione di alimenti.

**Prodotto raccomandato**



[Produttore di asciugamani per la pulizia viso](#)



[Asciugamano per il](#)



[Produttore di asciugamani per la pulizia](#)



[Telo mare in microfibra](#)

Se sei interessato ai nostri prodotti, non esitare a contattarci !!

Ricorda di guardare l'altro [asciugamano in microfibra](#) o [telo mare](#) Grazie! ^^